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





Assunto NC	Pré-tratamento	Pós-tratamento
mais profundo	NEWCASTLE CLINIC 26/10/19 10:22:48 PPJ NC120970, 12/09/70 MI 1.1 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR DEEPEST 1 L 11.35 mm	NEWCASTLE CLINIC 26/10/19 10:49:52 PPJ NC120970, 12/09/70 MI 1.2 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR POST TREATMENT 1 L 9.95 mm
longitudinal	NEWCASTLE CLINIC 26/10/19 10:24:45 PPJ NC120970, 12/09/70 MI 1.1 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR LONG 1 L 7.20 mm	NEWCASTLE CLINIC 26/10/19 10:52:47 PPJ NC120970, 12/09/70 MI 1.2 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR POST TREATMENT 1 L 4.65 mm
profundo	NEWCASTLE CLINIC 26/10/19 10:23:29 PPJ NC120970, 12/09/70 MI 1.1 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR DEEP 1 L 11.22 mm	NEWCASTLE CLINIC 26/10/19 10:52:01 PPJ NC120970, 12/09/70 MI 1.2 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR POST TREATMENT 1 L 5.13 mm
transversal	NEWCASTLE CLINIC 26/10/19 10:24:45 PPJ NC120970, 12/09/70 MI 1.1 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR TRANV 1 L 7.20 mm	NEWCASTLE CLINIC 26/10/19 10:52:47 PPJ NC120970, 12/09/70 MI 1.2 Tis 0.1 ML6-15 MSK Gen HYSTERECTOMY SCAR POST TREATMENT 1 L 4.65 mm
vascularidade	NEWCASTLE CLINIC 26/10/19 10:28:10 PPJ NC120970, 12/09/70 MI 1.0 Tis 0.8 ML6-15 MSK Gen HYSTERECTOMY SCAR VASC	NEWCASTLE CLINIC 26/10/19 10:55:11 PPJ NC120970, 12/09/70 MI 1.1 Tis 0.8 ML6-15 MSK Gen HYSTERECTOMY SCAR POST TREATMENT

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Alastair McLoughlin
www.McLoughlin-Scar-Release.com

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Re: D D DOB: 20.03.73

Ultrasound Examination – Anterior abdominal wall

Clinical Details: History of caesarean section scar in the lower abdominal wall.

Report: The anterior abdominal wall scar in the subcutaneous fat was barely visible and was difficult to measure. The approximate deepest dimension of the scar before treatment is 19.3mm but after the treatment decreased to 10.7mm.

The approximate depth of the scar which was measured just right of the midline (right lateral) was approximately 16.6mm before treatment but after treatment the approximate depth just right of midline decreased to approximately 3.6mm.

The approximate dimensions of the scar in longitudinal and transverse dimensions is 9.4mm x 9.9mm respectively before treatment but following treatment the scar tissue measures approximately 6.7mm x 7.3mm in maximum approximate longitudinal and transverse dimensions respectively.

There was no evidence of any vascularity noted in the scar or around the scar before treatment but following treatment, blood supply around the scar was noted especially in the subcutaneous adipose tissue but there was no evidence of any vascularity noted in the scar following treatment.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE

t: 0191 281 2636
f: 0191 281 2393

reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO10000008


Company Registration Number: 5843596

Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/SR

Scan Date: 26.10.19

30th October 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: P E DOB: 17.01.70

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Caesarean section noted.

Report: On ultrasound examination, scar tissue measures approximately 18.4mm x 14.9mm in maximum approximate longitudinal and transverse dimensions respectively before treatment but following treatment, there was a decrease in the dimensions of the scar tissue. The scar tissue measures approximately 8.8mm x 11.3mm in maximum longitudinal and transverse dimensions respectively.

Before treatment, there was no evidence of any vascularity noted in and around the scar but following treatment, there was vascularity noted around the scar in the anterior fascia covering the anterior aspect of the rectus abdominus muscle.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE
t: 0191 281 2636
f: 0191 281 2393
reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO1000008
Company Registration Number: 5843596
Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/SR

Scan Date: 26.10.19

30th October 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: C W DOB: 20.07.69

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Caesarean section noted.

Report: There is evidence of lower abdominal wall caesarean section scar. The deepest dimensions of the anterior abdominal wall scar in the region of the caesarean section measures approximately 19.5mm before treatment but following treatment, the deepest dimension of the scar decreased to approximately 14.2mm only. The approximate depth of the scar before treatment was 16.1mm especially to the right of the midline but following treatment, the approximate depth of the scar decreased to 10.6mm. Approximate dimensions of the scar are 15.1mm x 12.7mm in maximum longitudinal and transverse dimensions respectively before treatment but following treatment, the approximate dimensions of the scar are 8.4mm x 8.5mm in maximum longitudinal and transverse dimensions respectively. On power Doppler interrogation there was minimal vascularity noted around the scar, but no evidence of any vascularity in the scar tissue. Following treatment, there was increase in the vascularity around the scar tissue but again, no evidence of any abnormal vascularity noted in the scar tissue following treatment.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE

t: 0191 281 2636
f: 0191 281 2393

reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO10000008


Company Registration Number: 5843596

Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/SR

Scan Date: 26.10.19

30th October 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: A B DOB: 12.05.66

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Caesarean section noted.

Report: There is evidence of healed scar noted in the suprapubic region in the lower abdominal wall related to healed caesarean section scar.

Approximately deepest dimension of the scar before treatment is 17mm which decreased to 15.7mm following treatment. The depth of the scar just right of midline is approximately 10.9mm which decreased to 9mm following treatment.

Approximate dimensions of the scar are 14.2mm x 13.5mm and maximum longitudinal and transverse dimension respectively before treatment but following treatment, the approximate dimensions are 8.2mm x 11mm and maximum longitudinal and transverse dimensions respectively.

There was no evidence of any vascularity noted around the scar before treatment but following treatment, there was evidence of minimal vascularity noted in the scar and around the scar on power Doppler interrogation.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE
t: 0191 281 2636
f: 0191 281 2393
reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO10000008
Company Registration Number: 5843596
Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/SR

Scan Date: 26.10.19

30th October 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: K H DOB: 10.09.77

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Caesarean section noted.

Report: The deepest dimension of the scar in the midline is approximately 7mm before treatment but following treatment, the deepest dimension of the scar decreased to approximately 5.8mm.

The dimension of the scar especially in its maximum depth just right of midline is approximately 5.6mm before treatment but following treatment, this dimension decreased to approximately 3.6mm.

The approximate dimensions of the scar are 8.9mm x 11mm in maximum transverse and longitudinal dimensions respectively before treatment but following treatment, the approximate dimensions of the scar are 6.3mm x 10.7mm and maximum transverse and longitudinal dimensions respectively.

On power Doppler interrogation, there was no evidence of any vascularity noted in the scar or around the scar but following treatment, there was evidence of vascularity noted around the scar including mildly increased vascularity in the scar itself. Please note that this is a deep fascial scar and there was no evidence of any subcutaneous scar tissue especially in the subcutaneous fat on the ultrasound examination.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE
t: 0191 281 2636
f: 0191 281 2393
reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO10000008
Company Registration Number: 5843596
Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/LE

Scan Date: 26.10.19

5th November 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: J C D.O.B. 10.05.83

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Tummy Tuck scar noted.

Report: The deepest dimension of the scar is approximately 10 mm before treatment. Treatment and the deepest dimension decreased to approximately 7.8 mm. The approximate dimensions of the scar just right of midline is 7.7 mm in its maximum depth which decreased to approximately 7.1 mm following treatment.

The dimensions of the scar tissue is approximately 8 mm x 5.8 mm in maximum longitudinal and transverse dimensions respectively before treatment.

After treatment, the approximate dimensions of scar are 7.6 mm x 4 mm in maximum longitudinal and transverse dimensions respectively.

Before treatment, there was no evidence of vascularity in the scar tissue and there was no evidence of any vascularity noted around the scar tissue on power Doppler interrogation. Following treatment, there was increased vascularity noted in the scar and around the scar.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE

t: 0191 281 2636

f: 0191 281 2393

reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No: NO10000008

Company Registration Number: 5843596

Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

Ref: PPJR/SR

Scan Date: 26.10.19

30th October 2019

Mr Alastair McLoughlin

Germany

Dear Mr McLoughlin

Re: S S DOB: 13.03.86

Ultrasound Examination – Anterior abdominal wall

Clinical Details: Caesarean section noted.

Report: Before treatment, the depressed dimension of the scar tissue is approximately 9.3mm which decreased to approximately 7.5mm following treatment of the scar. The approximate dimensions of the scar just right of midline is 6.3mm before treatment but following treatment, the approximate dimensions of the scar is 4.4mm.

The approximate measurements of the scar is 5.3mm x 3.4mm in maximum longitudinal and transverse dimension respectively following treatment, the approximate dimensions of the scar of 4.8mm x 5.7mm in maximum longitudinal and transverse dimensions respectively.

On power Doppler interrogation, there was no evidence of any vascularity noted in and around the scar but following treatment, there was evidence of vascularity noted around the scar which was essentially noted just superficial and anterior to the scar in the subcutaneous soft tissues.

Yours sincerely

Dr. P P J Raju
Consultant Radiologist

Newcastle Clinic Limited
4 Towers Avenue
Jesmond
Newcastle Upon Tyne
NE2 3QE
t: 0191 281 2636
f: 0191 281 2393
reception@newcastleclinic.co.uk

Registered by the Care Quality Commission No. NO1000008
Company Registration Number: 5843596
Registered Office: Oakapple House, 1 John Charles Way, Leeds, LS12 6QA

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© Alastair McLoughlin